



Extracalm[®]

soothes itching and cools the skin

What is Extracalm[®]?

Extracalm[®] is a gel spray against itching.

What is Extracalm[®] used for?

Extracalm[®] is used for itchy, painful, red or irritated skin (not for open wounds) caused by:

- Insect bites
- First-degree burns
- Sunburn
- Dry skin
- Dermatological procedures (e.g. peelings, laser treatment)
- Friction (intact skin susceptible to developing bedsores)
- Allergic reactions caused by contact with certain substances

Extracalm[®] is safe and easy to apply, even on children and pregnant women.

How to use Extracalm[®]?

1. Shake the spray can.
2. Hold the spray can 10-15 cm from the skin.
3. Hold the spray can upright and spray the skin until it is completely covered.

Extracalm[®] is quickly absorbed by the skin. No need to rub it in.

Apply Extracalm[®] regularly, 2 to 6 times per day until the symptoms have disappeared.

Hold the spray can further away if Extracalm[®] feels too cold.

In the case of sunburn, apply Extracalm[®] frequently to minimize peeling and stay out of the sun.

Keep your eyes closed if Extracalm[®] is applied to the face.

If the skin treated with Extracalm[®] becomes wet (i.e. after swimming), dab dry and apply Extracalm[®] again.

How does Extracalm[®] work?

Extracalm[®] gel spray soothes itching, cools and restores the moisture balance of the skin.

Special precautions when using Extracalm[®]

- Do not apply Extracalm[®] to an open wound. If Extracalm[®] is nevertheless applied to an open wound, a stinging sensation may occur. Rinse the wound with clean water or a wound cleanser.
- Avoid contact with the eyes. If Extracalm[®] comes into contact with the eyes, rinse the eyes with clean water and consult a physician.
- Extracalm[®] does not contain a sunscreen agent. Use appropriate protection when sunbathing.
- Side-effects such as irritations or allergic reactions are rare. Do not use Extracalm[®] if there is a known allergy to parabens or any of the ingredients in Extracalm[®].
- If your general condition worsens (e.g. fever or headache) or if your skin condition worsens, consult a physician.
- Do not use Extracalm[®] for more than 30 days.
- Keep Extracalm[®] out of the reach of children.
- Keep away from sources of ignition - No smoking.
- Do not spray on a naked flame or any incandescent material.
- Extracalm[®] is a pressurized container. Protect from sunlight and do not expose to temperatures exceeding 50° C.
- Do not pierce or burn the spray can, even after use.
- Do not inhale the spray.
- Use the spray can only in well-ventilated areas. Insufficient ventilation can cause the formation of an explosive mix of fumes and air.
- Extracalm[®] is not harmful to the ozone layer.

Storage conditions

Keep at room temperature.

Do not use Extracalm[®] after the expiry date shown on the packaging. The expiry date (USE BY year/month) is the last day of the indicated month and year.

Composition

Purified water (Eur. Ph.), BCFA, Dimethicone, Acid colloidal polymer, Arginine, Methyl-p-hydroxybenzoate, Propyl-p-hydroxybenzoate, Macrogol.

Presentation*

- Spray can of 150 ml

**not available in every country*

Extracalm[®] is a medical device.